Device description

PuraStat is a synthetic haemostatic material in the form of a prefilled syringe, filled with a clear, 2.5% (range 1.8 to 3.0%) aqueous peptide solution, sterilized by aseptic filtration. The outer surface of the syringe and inner surface of the blister pack are sterilized by ethylene oxide.

Contact between the product and liquid such as blood causes the acidic peptide solution to be neutralized or become alkaline and, as a result, the peptide molecule, which has a β structure, quickly forms fibres in the aqueous solution, yielding a peptide hydrogel. The hydrogel quickly coats the point of bleeding and, by physically closing the superficial part of the broken blood vessel, causes blood coagulation in the deeper part of the vascular wall, resulting in haemostasis.

NOTE: Any product remaining in the body after haemostasis has been achieved, and the surgical procedure is complete, will be absorbed over time although some residue may remain for longer than 30 days.

Indications for use

PuraStat is indicated for haemostasis in the following situations encountered during surgery, when haemostasis by ligation or standard means is insufficient or impractical:

- Exudative haemorrhage from blood vessels and parenchyma of solid organs
- Exudative haemorrhage from vascular anastomoses
- Exudative haemorrhage from small vessels and capillary vessels of the GI tract

Circumstances and examples of these situations include the following:

Type of bleeding	Relevant clinical circumstances	Examples
Exudative haemorrhage from blood vessels and parenchyma of solid organ	Haemorrhagic point is a solid organ resection surface. Organ has a high-density vascular network and extra-vascular blood supply, e.g. liver or spleen. Any partial resection of such organs or removal of any organ attached to them.	Hepatectomy; splenectomy; cholecystectomy
Exudative haemorrhage from vascular anastomoses	Haemorrhagic point is a gap in the anastomoses. Blood with low clotting function (due to anticoagulation). Any vascular anastomosis, vessel to vessel or graft to vessel, or arteriotomy closure	Anastomosis to native or artificial vessel; surgery of the aorta or any peripheral arteries; coronary bypass; femoral bypass
Exudative haemorrhage from small vessels and capillary vessels of the GI tract	Haemorrhagic point is the mucosal resection site. Presence of digestive fluids (e.g. Bile, pancreatic fluid). Laparoscopic or endoscopic (transcatheter application).	Endoscopic mucosal resection of GI tract; Endoscopic submucosal dissection of upper GI tract; Laparoscopic resection of GI tract

Contra-indications

No contra-indications have been identified relating to the indications for use. The product has not been tested in pregnant women or in children.

Adverse events

The following problems and adverse events may occur when the product is used:

- Inflammatory responses, such as adhesions
- Other problems and adverse events, not mentioned here, may also occur.

Procedure for use

1. Before use

Read these instructions carefully.

PuraStat is made available in three fill volume variants. Select the appropriate variant relevant to the procedure to be completed.

Before using the product confirm that there is no damage to the package or the syringe. Stop using the product if any damage or leakage is noted. Connect the syringe to the optional nozzle included in the package, if appropriate to the type of surgery.

2. Method of use

- (A) Remove as much blood as possible from haemorrhagic site by for example applying adequate pressure
- (B) Apply an adequate quantity of PuraStat to the haemorrhagic site.
- (C) Repeat the application of PuraStat several times until haemostasis is achieved, removing any gelled PuraStat by gloved hand or with appropriate instruments between applications.
- (D) After haemostasis is confirmed, remove excessive gelled PuraStat from the application site by gloved hand or with instruments.

3. After use

If the optional dedicated nozzle has been used with the syringe, ensure that it remains connected on completion of the procedure.

Dispose of the syringe, nozzle and any unused PuraStat by placing into an appropriate container as clinical waste.

Warnings and precautions for use

1. PuraStat should be used only by surgeons and only if suturing, ligation or other standard means of haemostasis are insufficient or impractical. Note that the product is not intended to replace surgical haemostasis procedures such as suturing or ligation.
2. PuraStat is not intended for use on spurting and/or gushing haemorrhages. Such haemorrhages can usually be stopped by conventional means such as ligation and clipping.
3. Do not use PuraStat as the main haemostatic treatment in patients with a blood coagulation disorder.
4. Exposure of PuraStat to a saline solution or to blood before application to the haemorrhagic site decreases its haemostatic effect. Before applying PuraStat, remove saline solution and blood as much as possible from the haemorrhagic site.
5. Be careful not to contaminate PuraStat during use. Handle PuraStat aseptically.
6. Use the product promptly after opening the package to prevent contamination.
7. PuraStat is a sterile product. Do not re-sterilize PuraStat.
8. Dispose of PuraStat in an appropriate clinical waste container.
9. PuraStat is a single-use product and must not be reused. Each syringe and nozzle is intended to be used once only for a single patient to achieve haemostasis within a single surgical procedure only. Once expressed from the syringe and in contact with fluids such as blood, the PuraStat peptide solution will begin to change character and self-assemble thereby rendering it non-reusable. Any attempt to retain any unused portion remaining in the syringe after first use has the potential to transfer blood/tissue borne diseases and infection. Any attempt to resterilize the product will result in degradation and loss of performance.
10. PuraStat is not intended for injection.
11. Do not connect the PuraStat syringe to any device other than the dedicated nozzle supplied with the syringe. Do not attempt to connect the PuraStat syringe to any product that has not been specifically designed and indicated for use with the PuraStat syringe. When the dedicated nozzle is used with the syringe, ensure that nozzle does not become detached during procedure.
12. Note that PuraStat may cause discolouration of the surface of vascular grafts due to haemolysis of adherent erythrocytes which may appear as brown stains.

Storage

1. PuraStat should be stored in a refrigerator (from 2°C to 8°C).
2. Keep dry.

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